USER MANUAL

US/ES

11 GPM HOG SLAT INJECTOR

CONTENTS / CONTENIDO

SAFETY INSTRUCTIONS / CONSIGNAS DE SEGURIDAD

TECHNICAL SPECIFICATIONS/ CARACTERISTICAS

INSTALLATION GUIDE / GUÍA DE INSTALACIÓN

MAINTENANCE/ MANTENIMIENTO

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SAFETY INSTRUCTIONS / CONSIGNAS DE SEGURIDAD

Please read the following recommendations carefully before putting your dispenser into service. Then keep these safety instructions carefully in a location easily accessible to everyone who will need to use this unit.

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Opinion

This Manual cannot be reproduced, transmitted, transcribed and stored in a server or translated into any language, in full or in part, regardless of the means, without previous written agreement from the Manufacturer. The Manufacturer reserves the right to modify the characteristics of the product described in this user's manual at any time. The MANUFACTURER cannot be held responsible for damage resulting from use of your injector.

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Warning

You must comply with the safety instructions when installing, using and performing maintenance on the injector: use appropriate tools, protective clothing and safety goggles when working on the equipment and install the dosing unit so that it can be used without presenting a risk. Follow the instructions in this manual and take the appropriate safety measures for the water temperature and the type of liquid being drawn up. Take extreme care when using dangerous substances (solvents, acids, and corrosive, toxic, caustic and inflammable substances, etc.).

For dosing of these substances, consult with your distributor before use in order to confirm they are compatible with your injector.

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Max Pressure / Presión máxima: 85 PSI

Max Temperature / Temperatura máxima: 104°F

The injector should be installed in a position protected from frost.

The injector is delivered with: a mounting bracket/ a flexible suction tube with a strainer.

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
INSTALLATION GUIDE / GUÍA DE INSTALACIÓN

⚠ Warning

Personnel in charge of installing, using and maintaining this equipment must read this manual. When connecting an installation either to the public water supply or to its own water source, you must respect the regulations in force concerning protection of the source i.e. backflow prevention, etc. **IMPORTANT: The water must flow through the injector in the direction of the arrow molded on the injector body.**

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Installing the unit on a by-pass enables clean water to be supplied without operating the injector. Do not install your injector on the suction side of the supply pump (risk of siphoning).

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It is recommended that you should place an anti-siphon valve in installations in which there is a risk of siphoning.

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IMPORTANT ! It is the responsibility of the owner/operator to check that the flow and pressure of the installation do not exceed the injector characteristics.

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
It is the responsibility of the owner/operator to determine the correct amount of solution and injection ratio to obtain the desired result.

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The location of the injector and concentrate container should be accessible, but should never present a risk of pollution or contamination. It is recommended to label all water lines with a warning about the injected solution i.e. **IMPORTANT ! Not For Human Consumption**.

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The above flows may vary. This chart should only be used as a guide.

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MAINTENANCE / MANTENIMIENTO

Before any dismantling of the injector it should be rinsed to eliminate all traces of concentrated products and should not be left under pressure.

Never use metal tools. All the components are assembled manually. Never use excessive force. The injector should be rinsed regularly. We recommend that, each time the container of concentrate is changed, the suction tube and its strainer is plunged into a bucket of clean water and the unit is operated for 5 or 6 minutes.

After rinsing, transfer the suction tube to the new container and operate the unit until the tube is filled with concentrate. This rinsing procedure is imperative when changing the type or brand of product.

If the unit is not used over a long period, it should be rinsed and should not be left under pressure.

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
Turn off water supply. Disconnect unit from water source.

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To remove top seal from its seating, move the dosing piston up and down as far as it will go, as indicated by the arrows. This will dislodge the seal.

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To fit a new seal it is advisable to dismantle the valve, valve spring and pin. The seal must be positioned with the lip towards the top. When reassembling, be careful to respect the reverse order.

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Check that the pin is accurately centered and does not project at either end.

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- 3 Unscrew the black stem cap at the bottom of the injection stem. Change the seal. Refit the stem cap and tighten by hand without forcing. Raised center of seal towards the top.

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TROUBLESHOOTING GUIDE / MANUAL DE MANTENIMIENTO

This Injector was designed to give you trouble free and dependable service. However, if a minor problem occurs, we have prepared this troubleshooting guide for your use.

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Hog Slat, Inc., provides a limited warranty on the Hog Slat injector. The Term of Warranty is for a period of one-year from date of original purchase provided all warranty terms and conditions are followed.

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This Limited Warranty covers your unit for manufacturing and materials defects only. Warranty repair of the medicator will be done at the discretion of Hog Slat, Inc.

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There is no warranty expressed or implied relating in any way to products used in conjunction with the Hog Slat injector.

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There are no warranties, expressed or implied, which extend beyond those described above.

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Returns / Devoluciones

The return for warranty repair are made at the discretion of Hog Slat, Inc. Freight charges incurred when sending warranty repairs are the responsibility of the end-user. If the Warranty is valid, Hog Slat, Inc. pays freight costs when returning the repaired injector to end-users within the continental USA. If the warranty is not valid or if the end-user is outside the continental USA, freight is FOB. Do not ship a unit for repair until you have called one of the Hog Slat Service Centers and obtain a valid RA (Return Authorization) number.

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Include the following information when returning the unit: return address, billing address, name & telephone number of a contact person, reason for return and description of the problem. Write the RGA on the outside of the shipping box. Returned units MUST be drained of any chemicals, rinsed with clear water, and all water drained from the unit prior to packing it in the box for return.

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Return the unit only. Do not return hoses, filters, or fitting connected to the unit. Pack the unit well with plenty of padding to protect it from shipping damage.

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Any estimate of repair costs will be given upon request.

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The complete model reference and the serial number of your unit is stamped on the pump body.

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Please record this number in the space below and refer to it when you call your distributor for information, parts, and service.

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